



TROAP mouse mAb

Catalog No	YP-mAb-09960
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TROAP
Immunogen	Synthesized peptide derived from human TROAP AA range: 532-582
Specificity	This antibody detects endogenous levels of human TROAP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tastin;Trophinin assisting protein;Trophinin associated protein;Trophinin-assisting protein;Trophinin-associated protein;TROAP;Tastin (Trophinin-assisting protein) (Trophinin-associated protein)
Observed Band	110 kDa
Calculated Molecular Weight	778 aa, 84 kDa
Cell Pathway	Cytoplasm.
Tissue Specificity	Strong expression at implantation sites. Was exclusively localized to the apical side of the syncytiotrophoblast. Also found in macrophages.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images