



LENG1 mouse mAb

Catalog No	YP-mAb-09961
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LENG1
Immunogen	Synthesized peptide derived from human LENG1 AA range: 27-77
Specificity	This antibody detects endogenous levels of human LENG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LENG1;Leukocyte receptor cluster member 1
Observed Band	31 kDa
Calculated Molecular Weight	264 aa, 31 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images