



YRDC mouse mAb

Catalog No	YP-mAb-09962
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	YRDC DRIP3 IRIP
Immunogen	Synthesized peptide derived from human YRDC AA range: 214-264
Specificity	This antibody detects endogenous levels of human YRDC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DRIP3;hIRIP;IRIP;RP11 109P14.4;SUA5;YRDC;YRDC DRIP3 IRIP;YrdC domain-containing protein;mitochondrial (Dopamine receptor-interacting protein 3) (Ischemia/reperfusion-inducible protein homolog) (hIRIP); mitochondrial (Dopamine receptor-interacting protein 3) (Ischemia/reperfusion-inducible protein homolog) (hIRIP)
Observed Band	29 kDa
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Mitochondrion . Cell membrane ; Peripheral membrane protein . A large fraction localizes in the cytoplasm, whereas a smaller fraction is imported to mitochondria. .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images