



ADH1A mouse mAb

Catalog No	YP-mAb-09963
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ADH1A ADH1
Immunogen	Synthesized peptide derived from human ADH1A AA range: 38-88
Specificity	This antibody detects endogenous levels of human ADH1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.1.1.1;ADH1;ADH1A;ADH1A ADH1;Alcohol dehydrogenase 1A (EC 1.1.1.1) (Alcohol dehydrogenase subunit alpha)
Observed Band	40 kDa
Calculated Molecular Weight	375 aa, 40 kDa
Cell Pathway	Cytoplasm.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images