



FCF1 mouse mAb

Catalog No	YP-mAb-09965
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FCF1 C14orf111 CGI-35
Immunogen	Synthesized peptide derived from human FCF1 AA range: 56-106
Specificity	This antibody detects endogenous levels of human FCF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Bka;UTP24;CGI-35;C14orf111;FCF1;FCF1 C14orf111 CGI-35;rRNA-processing protein FCF1 homolog
Observed Band	23kDa
Calculated Molecular Weight	23kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images