



NAF1 mouse mAb

Catalog No	YP-mAb-09967
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NAF1
Immunogen	Synthesized peptide derived from human NAF1 AA range: 409-459
Specificity	This antibody detects endogenous levels of human NAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hNAF1;H/ACA ribonucleoprotein complex non-core subunit NAF1;NAF1;H/ACA ribonucleoprotein complex non-core subunit NAF1 (hNAF1)
Observed Band	70 kDa
Calculated Molecular Weight	54 kDa
Cell Pathway	Cytoplasm. Nucleus. Shuttles between the cytoplasm and the nucleus. Absent from the nucleolus (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images