



GOG6A mouse mAb

Catalog No	YP-mAb-09969
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GOLGA6A GLP GOLGA6
Immunogen	Synthesized peptide derived from human GOG6A AA range: 123-173
Specificity	This antibody detects endogenous levels of human GOG6A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GOG6A;GOLGA6A GLP GOLGA6;Golgin subfamily A member 6A (Golgin linked to PML) (Golgin-like protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Highly expressed in seminiferous tubes in testis. Highly expressed in spermatids, barely detectable in late pachytene spermatocytes, and not detectable in spermatogonia. Detected at intermediate levels in pancreas and lymph nodes, and at much lower levels in spleen, peripheral blood leukocytes, skeletal muscle, liver, lung, placenta, brain and heart.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images