



ST1C3 mouse mAb

Catalog No	YP-mAb-09970
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SULT1C3
Immunogen	Synthesized peptide derived from human ST1C3 AA range: 87-137
Specificity	This antibody detects endogenous levels of human ST1C3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ST1C3;SULT1C3;Sulfotransferase 1C3 (ST1C3) (EC 2.8.2.2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	[Isoform 1]: Not detectable in any of the tissues tested. ; [Isoform 2]: Expressed in the small intestine.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images