



C1AS1 mouse mAb

Catalog No	YP-mAb-09975
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C17orf76-AS1 C17orf45 NCRNA00188
Immunogen	Synthesized peptide derived from human C1AS1 AA range: 19-69
Specificity	This antibody detects endogenous levels of human C1AS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1AS1;C17orf76-AS1 C17orf45 NCRNA00188;Putative uncharacterized protein C17orf76-AS1;mitochondrial (C17orf76 antisense RNA 1) (C17orf76 antisense gene protein 1); mitochondrial (C17orf76 antisense RNA 1) (C17orf76 antisense gene protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images