



S15A2 mouse mAb

Catalog No	YP-mAb-09976
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC15A2 PEPT2
Immunogen	Synthesized peptide derived from human S15A2 AA range: 484-534
Specificity	This antibody detects endogenous levels of human S15A2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S15A2;SLC15A2 PEPT2;Solute carrier family 15 member 2 (Kidney H(+)/peptide cotransporter) (Oligopeptide transporter;kidney isoform) (Peptide transporter 2); kidney isoform) (Peptide transporter 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in kidney (PubMed:7756356). Not detected in intestine (PubMed:7756356).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images