



CNPD1 mouse mAb

Catalog No	YP-mAb-09977
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CNPPD1 C2orf24 CDABP0125 CGI-57
Immunogen	Synthesized peptide derived from human CNPD1 AA range: 207-257
Specificity	This antibody detects endogenous levels of human CNPD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNPD1;CNPPD1 C2orf24 CDABP0125 CGI-57;Protein CNPPD1 (Cyclin Pas1/PHO80 domain-containing protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images