



BABA1 mouse mAb

Catalog No	YP-mAb-09978
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BABAM1 C19orf62 MERIT40 NBA1 HSPC142
Immunogen	Synthesized peptide derived from human BABA1 AA range: 84-134
Specificity	This antibody detects endogenous levels of human BABA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MERIT40;BABAM1 C19orf62 MERIT40 NBA1 HSPC142;BRISC and BRCA1-A complex member 1 (Mediator of RAP80 interactions and targeting subunit of 40 kDa) (New component of the BRCA1-A complex)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Localizes at sites of DNA damage at double-strand breaks (DSBs). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images