



ZMYM1 mouse mAb

Catalog No	YP-mAb-09979
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZMYM1
Immunogen	Synthesized peptide derived from human ZMYM1 AA range: 904-954
Specificity	This antibody detects endogenous levels of human ZMYM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZMYM1;Zinc finger MYM-type protein 1
Observed Band	
Calculated Molecular Weight	126kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images