



DB127 mouse mAb

Catalog No	YP-mAb-09980
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DEFB127 C20orf73 DEFB27 UNQ1956/PRO6071
Immunogen	Synthesized peptide derived from human DB127 AA range: 44-94
Specificity	This antibody detects endogenous levels of human DB127
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DB127;DEFB127 C20orf73 DEFB27 UNQ1956/PRO6071;Beta-defensin 127 (Beta-defensin 27) (DEFB-27) (Defensin;beta 127); beta 127)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images