



CC062 mouse mAb

Catalog No	YP-mAb-09981
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C3orf62
Immunogen	Synthesized peptide derived from human CC062 AA range: 173-223
Specificity	This antibody detects endogenous levels of human CC062
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC062;C3orf62;Uncharacterized protein C3orf62
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Testis. Down-regulated in men with spermatocyte arrest.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images