



PRCD mouse mAb

Catalog No	YP-mAb-09982
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRCD
Immunogen	Synthesized peptide derived from human PRCD AA range: 4-54
Specificity	This antibody detects endogenous levels of human PRCD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRCD;Progressive rod-cone degeneration protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium, photoreceptor outer segment . Membrane ; Lipid-anchor ; Cytoplasmic side . Endoplasmic reticulum . Golgi apparatus . Localizes to photoreceptor disk membranes in the photoreceptor outer segment (PubMed:27613864). The secretion in media described in PubMed:24992209 is probably an experimental artifact (PubMed:24992209). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images