



UBIM mouse mAb

Catalog No	YP-mAb-09987
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAU
Immunogen	Synthesized peptide derived from human UBIM AA range: 18-68
Specificity	This antibody detects endogenous levels of human UBIM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	40S ribosomal protein S30;asr1;FAU;FAU encoded ubiquitin like protein;FAU1; FBR MuSV associated ubiquitously expressed;Finkel Biskis Reilly murine sarcoma virus (FBR MuSV) ubiquitously expressed (fox derived);Finkel Biskis Reilly murine sarcoma virus (FBR MuSV) ubiquitously expressed;Finkel Biskis Reilly murine sarcoma virus ubiquitously expressed;FLJ22986;Fub1;Fubi; MNSFbeta;Monoclonal nonspecific suppressor factor beta;Ribosomal protein S30;RPS30;S30;UBIM_HUMAN;Ubiquitin like protein fubi and ribosomal protein S30;Ubiquitin like protein fubi;Ubiquitin like S30 fusion protein;Ubiquitin-like protein FUBI;UBIM
Observed Band	10kD
Calculated Molecular Weight	8kD
Cell Pathway	
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images