



DGCR6 mouse mAb

Catalog No	YP-mAb-09988
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DGCR6
Immunogen	Synthesized peptide derived from human DGCR6 AA range: 34-84
Specificity	This antibody detects endogenous levels of human DGCR6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DGCRK6;C22orf12;D16H22S788E;D16Wis2;DGCR 8;Dgcr8;DGCR8 microprocessor complex subunit;DGCR8_HUMAN;DGCRK 6;DiGeorge syndrome critical region 8;DiGeorge syndrome critical region gene 8;Gy1; Microprocessor complex subunit DGCR8;pasha;DGCR6;Protein DGCR6 (DiGeorge syndrome critical region 6)
Observed Band	100kDa
Calculated Molecular Weight	86kD
Cell Pathway	Nucleus . Predominantly nuclear.
Tissue Specificity	Found in all tissues examined with highest expression in liver, heart and skeletal muscle. Lower levels in pancreas and placenta. Weak expression in brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images