



RDH8 mouse mAb

Catalog No	YP-mAb-09989
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RDH8 PRRDH
Immunogen	Synthesized peptide derived from human RDH8 AA range: 14-64
Specificity	This antibody detects endogenous levels of human RDH8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Photoreceptor outer segment all-trans retinol dehydrogenase;PRRDH;RDH8; RDH8_HUMAN;Retinol dehydrogenase 8 (all-trans);Retinol dehydrogenase 8; SDR28C2;Short chain dehydrogenase/reductase family 28C;member 22;RDH8 PRRDH;Retinol dehydrogenase 8 (EC 1.1.1.300) (Photoreceptor outer segment all-trans retinol dehydrogenase)
Observed Band	
Calculated Molecular Weight	34kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Detected in photoreceptor outer segments in the retina (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images