



TAF9B mouse mAb

Catalog No	YP-mAb-09990
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TAF9B TAF9L
Immunogen	Synthesized peptide derived from human TAF9B AA range: 153-203
Specificity	This antibody detects endogenous levels of human TAF9B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DN 7;DN7;TAF9B;TAF9L;TAFII31L;TFIID 31;TAF9B TAF9L;Transcription initiation factor TFIID subunit 9B (Neuronal cell death-related protein 7) (DN-7) (Transcription initiation factor TFIID subunit 9-like) (Transcription-associated factor TAFII31L)
Observed Band	30 kDa
Calculated Molecular Weight	28 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images