



ODAM mouse mAb

Catalog No	YP-mAb-09991
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ODAM APIN
Immunogen	Synthesized peptide derived from human ODA M AA range: 23-73
Specificity	This antibody detects endogenous levels of human ODA M
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	APIN;ODAM;ODAM APIN;Odontogenic ameloblast-associated protein (Apin)
Observed Band	
Calculated Molecular Weight	153 aa, 17 kDa
Cell Pathway	Secreted . Cytoplasm . Nucleus .
Tissue Specificity	Expressed in the junctional epithelium of healthy teeth. In periodontitis, absent in the pocket epithelium of the diseased periodontium but is detected in the gingival crevicular fluid.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images