



CF050 mouse mAb

Catalog No	YP-mAb-09992
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C6orf50 NAG19
Immunogen	Synthesized peptide derived from human CF050 AA range: 19-69
Specificity	This antibody detects endogenous levels of human CF050
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CF050;C6orf50 NAG19;Putative uncharacterized protein C6orf50 (Nasopharyngeal carcinoma-associated gene 19 protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images