



CRBA4 mouse mAb

Catalog No	YP-mAb-09993
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CRYBA4
Immunogen	Synthesized peptide derived from human CRBA4 AA range: 142-192
Specificity	This antibody detects endogenous levels of human CRBA4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Beta A4 crystallin;Beta crystallin A4;Beta-A4 crystallin;Beta-crystallin A4; CRBA4_HUMAN;CRYBA4;Crystallin beta A4;Crystallin;beta polypeptide A4; CTRCT23;Eye lens structural protein;MCOPCT4;CRBA4;Beta-crystallin A4 (Beta-A4 crystallin)
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images