



PIGB mouse mAb

Catalog No	YP-mAb-09994
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PIGB
Immunogen	Synthesized peptide derived from human PIGB AA range: 320-370
Specificity	This antibody detects endogenous levels of human PIGB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPI mannosyltransferase 3;GPI mannosyltransferase III;GPI MT III;PIG B;PIGB; GPI mannosyltransferase 3 (EC 2.4.1.-) (GPI mannosyltransferase III) (GPI-MT-III) (Phosphatidylinositol-glycan biosynthesis class B protein) (PIG-B)
Observed Band	50-55 kDa
Calculated Molecular Weight	65 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images