



VAT1L mouse mAb

Catalog No	YP-mAb-09995
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VAT1L KIAA1576
Immunogen	Synthesized peptide derived from human VAT1L AA range: 118-168
Specificity	This antibody detects endogenous levels of human VAT1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VAT1L;VAT1L KIAA1576;Synaptic vesicle membrane protein VAT-1 homolog-like (EC 1.-.-.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Detected in skin fibroblasts.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images