



MS4A5 mouse mAb

Catalog No	YP-mAb-09996
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MS4A5 CD20L2 TETM4
Immunogen	Synthesized peptide derived from human MS4A5 AA range: 46-96
Specificity	This antibody detects endogenous levels of human MS4A5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MS4A5;MS4A5 CD20L2 TETM4;Membrane-spanning 4-domains subfamily A member 5 (CD20 antigen-like 2) (Testis-expressed transmembrane protein 4)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Expressed at high level in the testis. Detected also in the pancreas, heart and in the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images