



ZSC21 mouse mAb

Catalog No	YP-mAb-10000
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZSCAN21 ZFP38 ZNF38
Immunogen	Synthesized peptide derived from human ZSC21 AA range: 222-272
Specificity	This antibody detects endogenous levels of human ZSC21
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZSC21;ZSCAN21 ZFP38 ZNF38;Zinc finger and SCAN domain-containing protein 21 (Renal carcinoma antigen NY-REN-21) (Zinc finger protein 38 homolog) (Zfp-38)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images