



NAT8 mouse mAb

Catalog No	YP-mAb-10001
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NAT8 CML1 GLA TSC501
Immunogen	Synthesized peptide derived from human NAT8 AA range: 38-88
Specificity	This antibody detects endogenous levels of human NAT8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Acetyltransferase 2;ATase2;Camello-like protein 1;CCNAT;CML1;NAT8;NAT8 CML1 GLA TSC501;Probable N-acetyltransferase 8 (EC 2.3.1.-) (Camello-like protein 1)
Observed Band	25 kDa
Calculated Molecular Weight	26 kDa
Cell Pathway	Endoplasmic reticulum-Golgi intermediate compartment membrane; Single-pass type II membrane protein. Endoplasmic reticulum membrane; Single-pass type II membrane protein.
Tissue Specificity	Preferentially expressed in liver and kidney. Also detected in brain (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images