



SPDL1 mouse mAb

Catalog No	YP-mAb-10005
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Synthesized peptide derived from human SPDL1 AA range: 317-367
Specificity	This antibody detects endogenous levels of human SPDL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC99;hSpindly;MU-RMS-40.4A;Protein Spindly;SPDL1;spdly;Putative speedy protein-like protein LOC442572
Observed Band	70 kDa
Calculated Molecular Weight	605 aa, 70 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images