



K0494 mouse mAb

Catalog No	YP-mAb-10006
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIAA0494
Immunogen	Synthesized peptide derived from human K0494 AA range: 271-321
Specificity	This antibody detects endogenous levels of human K0494
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F hand domain containing protein KIAA0494;EF-hand domain-containing protein KIAA0494;K0494_HUMAN;Kiaa0494;K0494
Observed Band	
Calculated Molecular Weight	55kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images