



MA1A2 mouse mAb

Catalog No	YP-mAb-10007
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MAN1A2 MAN1B
Immunogen	Synthesized peptide derived from human MA1A2 AA range: 508-558
Specificity	This antibody detects endogenous levels of human MA1A2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MA1A2;MAN1A2 MAN1B;Mannosyl-oligosaccharide 1;2-alpha-mannosidase IB (EC 3.2.1.113) (Mannosidase alpha class 1A member 2) (Processing alpha-1;2-mannosidase IB) (Alpha-1;2-mannosidase IB)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane; Single-pass type II membrane protein.
Tissue Specificity	Highest levels of expression in placenta and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images