



P4R3B mouse mAb

Catalog No	YP-mAb-10008
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SMEK2 KIAA1387 PP4R3B PPP4R3B
Immunogen	Synthesized peptide derived from human P4R3B AA range: 425-475
Specificity	This antibody detects endogenous levels of human P4R3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PSY2;smk1;FLFL2;SMEK2;PP4R3B;PPP4R3B;P4R3B;SMEK2 KIAA1387 PP4R3B PPP4R3B;Serine/threonine-protein phosphatase 4 regulatory subunit 3B (SMEK homolog 2)
Observed Band	97kDa
Calculated Molecular Weight	97kDa
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Nucleus . In interphase localized in the cytoplasm and (with higher levels) the nucleus. During metaphase located in pericentriolar regions.
Tissue Specificity	Moderately expressed in tissues and specific brain regions examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images