



TUTLB mouse mAb

Catalog No	YP-mAb-10014
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IGSF9B KIAA1030
Immunogen	Synthesized peptide derived from human TUTLB AA range: 910-960
Specificity	This antibody detects endogenous levels of human TUTLB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TUTLB;IGSF9B KIAA1030;Protein turtle homolog B (Immunoglobulin superfamily member 9B) (IgSF9B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane ; Single-pass type I membrane protein . Cell junction, synapse, postsynaptic density .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images