



COA1 mouse mAb

Catalog No	YP-mAb-10015
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COA1 C7orf44
Immunogen	Synthesized peptide derived from human COA1 AA range: 16-66
Specificity	This antibody detects endogenous levels of human COA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C7orf44;MITRAC15;Mitochondrial translation regulation assembly intermediate of cytochrome c oxidase protein of 15 kDa;Cytochrome c oxidase assembly factor 1 homolog;COA1;COA1 C7orf44;Cytochrome oxidase assembly protein 1 homolog
Observed Band	
Calculated Molecular Weight	146 aa, 17 kDa
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images