



VOPP1 mouse mAb

Catalog No	YP-mAb-10017
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VOPP1 ECOP GASP
Immunogen	Synthesized peptide derived from human VOPP1 AA range: 120-170
Specificity	This antibody detects endogenous levels of human VOPP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECOP;GASP;VOPP1;VOPP1 ECOP GASP;Vesicular;overexpressed in cancer; prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N); overexpressed in cancer; prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N)
Observed Band	
Calculated Molecular Weight	172 aa, 19 kDa
Cell Pathway	Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . When overexpressed, localizes in the nucleus and perinuclear regions.
Tissue Specificity	Widely expressed with highest levels in thymus and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images