



BIEA mouse mAb

Catalog No	YP-mAb-10020
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BLVRA BLVR BVR
Immunogen	Synthesized peptide derived from human BIEA AA range: 84-134
Specificity	This antibody detects endogenous levels of human BIEA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BIEA_HUMAN;Biliverdin IX alpha reductase;Biliverdin reductase A;Biliverdin-IX alpha-reductase;BLVR A;BLVR;Blvra;BVR A;BVR;BVRA;Zinc metalloprotein; BIEA;BLVRA BLVR BVR;Biliverdin reductase A (BVR A) (EC 1.3.1.24) (Biliverdin-IX alpha-reductase)
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images