



CQ085 mouse mAb

Catalog No	YP-mAb-10022
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C17orf85
Immunogen	Synthesized peptide derived from human CQ085 AA range: 253-303
Specificity	This antibody detects endogenous levels of human CQ085
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CQ085;C17orf85;Uncharacterized protein C17orf85 (Protein ELG)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images