



ST1E1 mouse mAb

Catalog No	YP-mAb-10025
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SULT1E1 STE
Immunogen	Synthesized peptide derived from human ST1E1 AA range: 209-259
Specificity	This antibody detects endogenous levels of human ST1E1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EST;EST 1;Estrogen sulfotransferase;ST1E1;STE;Sulfotransferase 1E1;SULT1E1;SULT1E1 STE;Estrogen sulfotransferase (EC 2.8.2.4) (EST-1) (Sulfotransferase 1E1) (ST1E1) (Sulfotransferase;estrogen-preferring); estrogen-preferring)
Observed Band	32-35 kDa
Calculated Molecular Weight	294 aa, 35 kDa
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Liver, intestine and at lower level in the kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images