



SE1L3 mouse mAb

Catalog No	YP-mAb-10026
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SEL1L3 KIAA0746
Immunogen	Synthesized peptide derived from human SE1L3 AA range: 857-907
Specificity	This antibody detects endogenous levels of human SE1L3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SE1L3;SEL1L3 KIAA0746;Protein sel-1 homolog 3 (Suppressor of lin-12-like protein 3) (Sel-1L3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images