



TIM16 mouse mAb

Catalog No	YP-mAb-10027
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PAM16 MAGMAS TIM16 TIMM16 CGI-136
Immunogen	Synthesized peptide derived from human TIM16 AA range: 46-96
Specificity	This antibody detects endogenous levels of human TIM16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E3 ubiquitin-protein ligase TRIM16;EBBP;EC:2.3.2.27;Estrogen-responsive B box protein;tripartite motif containing 16;TIM16;PAM16 MAGMAS TIM16 TIMM16 CGI-136;Mitochondrial import inner membrane translocase subunit TIM16 (Mitochondria-associated granulocyte macrophage CSF-signaling molecule) (Presequence translocated-associated motor subunit PAM16)
Observed Band	64 kDa
Calculated Molecular Weight	64 kDa
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images