



ZN746 mouse mAb

Catalog No	YP-mAb-10029
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF746 PARIS
Immunogen	Synthesized peptide derived from human ZN746 AA range: 37-87
Specificity	This antibody detects endogenous levels of human ZN746
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PARIS;Parkin-interacting substrate;Zinc finger protein 746;ZN746_HUMAN;ZNF746;ZN746;ZNF746 PARIS;Zinc finger protein 746 (Parkin-interacting substrate) (PARIS)
Observed Band	69kD
Calculated Molecular Weight	69kD
Cell Pathway	Cytoplasm . Nucleus . Mainly localizes to the cytoplasm; probably translocates to the nucleus to repress selected genes.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images