



TPC10 mouse mAb

Catalog No	YP-mAb-10032
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRAPPC10 EHOC1 TMEM1
Immunogen	Synthesized peptide derived from human TPC10 AA range: 468-518
Specificity	This antibody detects endogenous levels of human TPC10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPC10;TRAPPC10 EHOC1 TMEM1;Trafficking protein particle complex subunit 10 (Epilepsy holoprosencephaly candidate 1 protein) (EHOC-1) (Protein GT334) (Trafficking protein particle complex subunit TMEM1) (Transport protein partic
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, cis-Golgi network .
Tissue Specificity	Expressed in all tissues examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images