



GLOD4 mouse mAb

Catalog No	YP-mAb-10035
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GLOD4 C17orf25 CGI-150 My027
Immunogen	Synthesized peptide derived from human GLOD4 AA range: 53-103
Specificity	This antibody detects endogenous levels of human GLOD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C17orf25;CGI 150;CGI150;CGI-150;GLOD 4;GLOD4;GLOD4 C17orf25 CGI-150 My027;Glyoxalase domain-containing protein 4
Observed Band 34kD	35 kDa
Calculated Molecular Weight	289 aa, 33 kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	Expressed in heart, brain, liver, kidney, pancreas and placenta. Not expressed in skeletal muscle and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images