



MSL2 mouse mAb

Catalog No	YP-mAb-10036
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MSL2 KIAA1585 MSL2L1 RNF184
Immunogen	Synthesized peptide derived from human MSL2 AA range: 9-59
Specificity	This antibody detects endogenous levels of human MSL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MS L2;MSL 2;MSL2 like 1;MSL2L1;RNF184;MSL2;MSL2 KIAA1585 MSL2L1 RNF184;Male-specific lethal 2 homolog (MSL-2) (Male-specific lethal 2-like 1) (MSL2-like 1) (Male-specific lethal-2 homolog 1) (RING finger protein 184)
Observed Band	70-75 kDa
Calculated Molecular Weight	577 aa, 63 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images