



VAC14 mouse mAb

Catalog No	YP-mAb-10037
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VAC14 TAX1BP2 TRX
Immunogen	Synthesized peptide derived from human VAC14 AA range: 569-619
Specificity	This antibody detects endogenous levels of human VAC14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ArPIKfyve;Protein VAC14 homolog;Tax1 binding protein 2;TAX1BP2;TRX;VAC14;Vac14 homolog (S. cerevisiae);VAC14 TAX1BP2 TRX;Protein VAC14 homolog (Tax1-binding protein 2)
Observed Band	88 kDa, 65 kDa
Calculated Molecular Weight	88 kDa
Cell Pathway	Endosome membrane . Microsome membrane . Mainly associated with membranes of the late endocytic pathway.
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images