



THOC6 mouse mAb

Catalog No	YP-mAb-10039
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	THOC6 WDR58 PSEC0006
Immunogen	Synthesized peptide derived from human THOC6 AA range: 79-129
Specificity	This antibody detects endogenous levels of human THOC6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	fSAP35;THO complex subunit 6 homolog;THOC6;WDR58;THOC6 WDR58 PSEC0006;THO complex subunit 6 homolog (Functional spliceosome-associated protein 35) (fSAP35) (WD repeat-containing protein 58)
Observed Band	35 kDa
Calculated Molecular Weight	38 kDa
Cell Pathway	Nucleus . Nucleus speckle .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images