



B3GN5 mouse mAb

Catalog No	YP-mAb-10040
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	B3GNT5
Immunogen	Synthesized peptide derived from human B3GN5 AA range: 167-217
Specificity	This antibody detects endogenous levels of human B3GN5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B3GN5;B3GNT5;Lactosylceramide 1;3-N-acetyl-beta-D-glucosaminyltransferase (EC 2.4.1.206) (Lactotriaosylceramide synthase) (Lc(3)Cer synthase) (Lc3 synthase) (UDP-GlcNAc:beta-Gal beta-1;3-N-acetylglucosaminyltransf
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	Widely expressed. Highly expressed in lung, colon, placenta, testis, pituitary gland and cerebellum. Weakly expressed in brain, liver, spleen, lymph node and thymus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images