



GLT13 mouse mAb

Catalog No	YP-mAb-10041
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GALNT13 KIAA1918
Immunogen	Synthesized peptide derived from human GLT13 AA range: 485-535
Specificity	This antibody detects endogenous levels of human GLT13
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GLT13;GALNT13 KIAA1918;Polypeptide N-acetylgalactosaminyltransferase 13 (EC 2.4.1.41) (Polypeptide GalNAc transferase 13) (GalNAc-T13) (pp-GaNTase 13) (Protein-UDP acetylgalactosaminyltransferase 13) (UDP-GalNAc:polypeptide
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	Specifically expressed in neuronal cells. Expressed in fetal brain, whole adult brain, cerebral cortex and cerebellum. Not expressed in other tissues tested.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images