



LY66C mouse mAb

Catalog No	YP-mAb-10044
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LY6G6C C6orf24 G6C NG24 UNQ1947/PRO4430
Immunogen	Synthesized peptide derived from human LY66C AA range: 53-103
Specificity	This antibody detects endogenous levels of human LY66C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LY66C;LY6G6C C6orf24 G6C NG24 UNQ1947/PRO4430;Lymphocyte antigen 6 complex locus protein G6c
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor .
Tissue Specificity	Highly expressed at the leading edges of cells, on filopodia.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images