



PGAP3 mouse mAb

Catalog No	YP-mAb-10045
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PGAP3 CAB2 PERLD1 UNQ546/PRO1100
Immunogen	Synthesized peptide derived from human PGAP3 AA range: 252-302
Specificity	This antibody detects endogenous levels of human PGAP3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGLA546;CAB2;COS16 homolog;hCOS16;per1 like domain containing 1; PERLD1;PGAP3;PP1498;PGAP3 CAB2 PERLD1 UNQ546/PRO1100;Post-GPI attachment to proteins factor 3 (COS16 homolog) (hCOS16) (Gene coamplified with ERBB2 protein) (PER1-like domain-containing protein 1)
Observed Band	36-46 kDa
Calculated Molecular Weight	320 aa, 36 kDa
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Mainly localizes to Golgi apparatus. .
Tissue Specificity	Ubiquitously expressed, with highest levels in thyroid and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images